

Human Peroxisome Proliferator-Activated Receptor γ Transactivation Assay (hPPAR γ -TA)

References

Garoché C, Boulahtouf A, Grimaldi M, Toporova L, den Broeder MJ, Legler J, Bourguet W, Balaguer P. "Interspecies differences in activation of Peroxisome proliferator-activated receptor γ by pharmaceutical and environmental chemicals". *Environ Sci Technol* 2021. DOI: 10.1021/acs.est.1c04318.

Léger T, Balaguer P, Le Hégarat L, Fessard V. Fate and PPAR γ and STATs-driven effects of the mitochondrial complex I inhibitor tebufenpyrad in liver cells revealed with multi-omics. *J Hazard Mater*. 2023 Jan 15;442:130083. DOI: 10.1016/j.jhazmat.2022.130083.

Summary of the method

The peroxisome proliferator-activated receptor γ (PPAR γ) is a well-known regulator of adipocyte differentiation, immune responses, proliferation of different tissues, and metabolism including fat storage. PPAR γ has been the target of drug discovery projects for a long period. These medicines include glitazones which enhance insulin sensitivity. There are numerous other metabolites and xenobiotics able to interfere with binding to PPARs. However, due to the obligatory heterodimerization of PPAR with the retinoid X receptor (RXR), it is complicated to distinguish effects caused by PPAR ligands from effects caused by RXR ligands. Moreover, dysregulation of PPAR γ signaling is associated with the development of several tumors.

This in vitro method measures the specific activation or inhibition of chemical substances on the human PPAR γ through a luciferase reporter gene. Similar to other transactivation assays, PPARTA relies on the binding and activation of chemicals to PPAR γ to induce the luciferase. This can be turned into an antagonistic assay using rosiglitazone that may be competed by test chemicals.

Readout and Result

This test elucidates the Molecular Initiating Event of the activation/inhibition of hPPAR γ . The assay is based on HeLa cells expressing the DNA binding domain (DBD) of the yeast transcription factor GAL4 fused to the ligand binding domain (LBD) of human hPPAR γ along with the responsive element of GAL4 fused to luciferase was developed. This allows for transactivation assessment of activation of hPPAR γ measuring the luciferase activation, enabling identification of PPAR γ agonists and antagonists.




Figure 1: Workflow of the hPARTA. After seeding of the cells in 96-well plates, they are treated with a test compound for 24 hours. Thereafter, D-Luciferin is added, and Luminescence is measured.

Necessary equipment

- standard cell culture equipment (e.g. incubator, laminar flow hood)
- 96-wells white opaque culture plates (e.g. Greiner bio-one 655083-905, CellStar; Dutscher, Brumath, France)
- Cell cytotoxicity measurement (e.g. Alamar blue, or Neutral red).
 - DMEM/F-12 without phenol red (Gibco 21041-025) Test medium
 - 5% DCC-treated FBS
 - 1% v/v penicillin/streptomycin (Gibco 15070-63)

Culture medium	<ul style="list-style-type: none"> - DMEM/F-12 with phenol red (Gibco 31331-028) - 10% FBS (Eurobio CVFSVF00) - 1% v/v penicillin/streptomycin (Gibco 15070-63) - 1 mg/mL geneticin (Invivogen ant-gn)
Test medium	<ul style="list-style-type: none"> - DMEM/F-12 without phenol red (Gibco 21041-025) - 5% DCC-treated FBS - 1% v/v penicillin/streptomycin (Gibco 15070-63)
Luminescence medium	<ul style="list-style-type: none"> - DMEM/F-12 without phenol red (Gibco 21041-025) - 5% DCC-treated FBS - 1% v/v penicillin/streptomycin (Gibco 15070-63) - 0.3 mM D-luciferin (Perkin Elmer 122799)

- Luminometer
 - Example : MicroBeta Wallac luminometer (Perkin-Elmer). Reading of microplates, from above.
 - Reading is optimal 20 minutes after the addition of the luminescence medium.